

## OSSGEN Dualpor Collagen PUTTY **ANCHOR TYPE**



## Main Features

- Porous calcium phosphate bone graft material with 'Inter-linked Macro/Micro porous structure'
- 60% Hydroxyapatite / 40%  $\beta$ -Tricalcium phosphate. (BCP: Biphasic Calcium phosphate)
- Trabecular structure similar to human cancellous bone and stable bone resorption in the body.
- Blood vessel growth and bone cell migration are possible through the connected macro porous.
- Very high autogenous bone occupancy (about 80% processing rate  $\rightarrow$  80% autogenous bone occupancy).
- The micro-porous structure on the skeletal surface is very advantageous for bone cell attachment.
- The optimal Scaffold for tissue engineering bone regeneration.
- Excellent biocompatibility.
- Excellent cell formation and adhesion.
- Playback effect.
- Convenience of use. (Free molding)
- Hemostatic effect and anti-adhesion effect.
- Superior bone conductivity and higher autogenous bone occupancy than conventional bone graft materials. (up to 85%)

## Raw material & Structure



Safe Raw Material

OssGen

Porous Membrane Structure




- Produced in New Zealand where is BSE-free nation
- Extraction from Bovine Amniotic membrane which is under 24 months
- Quality inspection by thorough Quality Assurance



- Production method through ingenious own technology
- White-dense sponge structure
- Easy to store & release drug by porosity membrane




## Structural advantage



Specimen  
( $\mu$ CT)

Result for 3month specimen-Micro Computed Tomography



| S University Hospital Research Results |

Fig. 1

- A** The  $\mu$ CT image of specimen with only defect and bone defects is indicated by red box.
- B** A  $\mu$ CT image of specimen grated DBM on defect and it is indicated by red box.
- C** A  $\mu$ CT image of specimen grafted collagen putty on defect and it is indicated by red box.

The control group showed a lot of healing in the cortical bone, but there were empty spaces in the trabecular bone. However, the treat group is more healed in cortical and trabecular bone healing than the others.




Fig. 2

The following are the Bone Mineral Densito metry(BMD, mg/cc) and bone volume for each sample.

**Negative control** 1849.12, 104.74

**DBM group** 2021.04, 149.021


**Treat(grafted Collagen putty)** 1932.03, 96.1124

	Bone volume(mm^3)	BMD(mg/cc)
N.C	104.74	1849.12
DBM	96.1124	1932.03
Treat	149.021	2021.04

# Dualpor Collagen PUTTY Anchor Type **Character**


Interconnected  
Macro/  
Micro pore  
system

OSG  
DualPor



Osteoblast adhesion on Micro pore surface after  
osteoblast moves to Simultaneous bone  
regeneration in pore inside&outside

Commercials  
Products




Cell adhesion on bone  
graft material outside

Comparison  
of  
Commercial  
Products

Cell culture  
behavior


A  
Company  
(US)

- Osteoblast adhesion only on bone graft material outside by closed pore
- No osteoblast in closed pore inside




B  
Company  
(US)

- Osteoblast adhesion only on bone graft material outside by closed pore
- No osteoblast in closed pore inside




OssGen  
(KR)

- Homogenous osteoblast adhesion along endoskeleton by 'inter-linked' open pore structure
- New bone occupancy ratio : **80%**




In-vivo  
Histology

Beagle dog  
model  
(6 month)




- Old bone (Upper left)
- Bone graft material (Brown color)
- New bone (Blue color)
- Osteoid (Ivory color around bone graft material)
- **80% porosity ratio of bone graft material → 80% autogenous bone occupancy ratio (Higer among domestic and foreign products)**
- Simultaneous & homogeneous bone regeneration of bone graft materail inside & outside through blood vessel growth & osteoblast movement to the inside of Macro/Micro porosity

## Collagen Putty Type with Tissue Engineering Courtesy of Prof. E. Park, KNU



Immunohistochemical analysis of new bone at 10 weeks after surgery  
(New bone formation measured by H&E staining)  
(a, b) defect only (c, d) BCP only (e, f) BCP + ATSC  
(g, h) BCP + BMSC (i, h) BCP + ATSC + VEGF  
(k, l) BCP + BMSC + VEGF \*NB : New bone






Immunohistochemical analysis of new bone at 10 weeks after surgery  
(New bone formation measured by masson's Trichrome staining)  
(a, b) defect only (c, d) BCP only (e, f) BCP + ATSC  
(g, h) BCP + BMSC (i, h) BCP + ATSC + VEGF  
(k, l) BCP + BMSC + VEGF \*collagen matrix formation



# OSSGEN


## Total Solution

Excellent osseointegration  
between bone graft material and  
New Bone (Osseointegration)



| Before blood mxing | | After blood mxing |

### Structure & Size



Model Name	A Top mm	B Bottom mm	C Height mm	D Top hole mm	E Bottom hole mm	F Hole Height mm	Weight g
Dualpor Collagen Putty H-A	7.2	4	45	3	1	25	0.55
Dualpor Collagen Putty H-B	7.2	4	50	3	1	40	0.59
Dualpor Collagen Putty H-C	8.5	5	45	6	1	25	0.72
Dualpor Collagen Putty H-D	8.5	5	50	6	1	40	0.73
Dualpor Collagen Putty H-E	6	5	40	3	1	25	0.46
Dualpor Collagen Putty H-F	6	5	50	3	1	30	0.57
Dualpor Collagen Putty H-G	7	5.5	40	3.3	1	25	0.59
Dualpor Collagen Putty H-H	7	5.5	50	3.3	1	30	0.75
Dualpor Collagen Putty H-I	6	5	40	3	1	30	0.45
Dualpor Collagen Putty H-J	6	5	50	3	1	30	0.57
Dualpor Collagen Putty H-K	7	5.5	40	3.3	1	30	0.58
Dualpor Collagen Putty H-L	7	5.5	50	3.3	1	30	0.75
Dualpor Collagen Putty H-M	7.2	4	45	3	1	30	0.55
Dualpor Collagen Putty H-N	7.2	4	50	3	1	30	0.61
Dualpor Collagen Putty H-O	8.5	5	45	6	1	30	0.69
Dualpor Collagen Putty H-P	8.5	5	50	6	1	30	0.78
Dualpor Collagen Putty H-Q	6	5	40	3	1	39	0.43
Dualpor Collagen Putty H-R	6	5	50	3	1	49	0.54
Dualpor Collagen Putty H-S	7	5.5	40	3.3	1	39	0.57
Dualpor Collagen Putty H-T	7	5.5	50	3.3	1	49	0.71
Dualpor Collagen Putty H-U	7.2	4	45	3	1	44	0.52
Dualpor Collagen Putty H-V	7.2	4	50	3	1	49	0.58
Dualpor Collagen Putty H-W	8.5	5	45	6	1	44	0.67
Dualpor Collagen Putty H-X	8.5	5	50	6	1	49	0.61



423-13, Hyeoksin-daero, Dong-gu, Daegu, Republic of Korea

Tel. +82-53-811-8191

Fax. +82-53-811-8192

E-mail. [daniel@ossenglobal.co.kr](mailto:daniel@ossenglobal.co.kr)